INTERNACIONAL FARMACEUTICA, S.A. DE C.V.

INTERNACIONAL FARMACEUTICA, S.A. DE C.V.


Our company fully complies with the
ISO 9001:2000 and ISO 13485:2003,
International Quality Standards, as well as
the Good Manufacturing Practice (GMP) guidelines.
Moreover, our 40 years experience in the medical devices
field allows us to insure premium quality
for all our products.

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



Depending on the sharpening technique employed, our needles can be classified as follows:

Standard or Normal sharpened Needles, PREMIUM Microsharpened Needles, ATRALLOY Needles, and our Supra-sharp Needles.

Our needles can be described according to the following characteristics:

- Curvature
- Length
- Wire Diameter
- Type of needle tips

In the pages that follow you will find our extensive catalogue of needles. Any specific needle which is not illustrated can be individually requested.

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


300 STAINLESS STEEL plus

300 STAINLESS STEEL plus

Atramat® has proudly developed a new generation of needles which perfectly combine a special stainless steel series alloy, with optimum geometric design, and a technologically innovative coating, resulting in an ideal harmonic needle.

Alloy

The technological development adopted for our new line of precision surgical needles is subject to the optimal selection of steel alloy based on stringent criteria. This is to ensure that the end product is of superior quality and suitably designed for the most complex specialty surgery.

Geometry

The exceptional geometric design of both the body and tip of the needle, means that it offers surgeons both outstanding strength while penetrating tissue and precise control.

Coating

A special thermal coating results in superior and easier penetration of tissue, coupled with minimal resistance and smoothness. This in turn reduces the effect of dragging, and means that the needle point remains sharper for a longer duration.

- Proven effectiveness.
- Technologically innovative.

ATRALLOY®, NON-REFLECTIVE DARK NEEDLE


Besides of the remarkable characteristics of the ATRALLOY® Needle in its alloy, geometry and coating, we have added a dark non-reflective finish for an optimum contrast in the surgical field, avoiding light reflection and enhancing the visualization. This is very convenient in micro-vascular surgery.

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This tip profile is designed to provide easy penetration of intended tissue. It is widely used in tissue with moderate ease of penetration, such as the peritoneum-fascia muscle, fatty tissue and intestine. It is also widely used in veins and arteries. Specialities include: General Surgery, Gastro-intestinal, Cardiovascular and Gynaecology.

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
The tapered cross section immediately behind the tip has been flattened to an oval shape in our round bodied needles; facilitating the penetration through tissue layers. It is widely used in General and Gastro-intestinal Surgery.

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
Cutting needle characterised by its triangular cutting pattern, with the apex cutting edge on the outside of the needle curvature. Widely used in Reconstructive Plastic Surgery and cuticular closing.

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
Cutting needle characterised by its triangular cutting pattern, with the apex cutting edge on the inside of the needle. Widely used in Reconstructive Plastic Surgery and in cuticular closing.

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
Widely used for tissue with tough penetration that requires a circular pattern. This avoids damage to the tissue from the tension caused by the suture. It is also used for tendons, Orthopaedics and sternum closure, due to its great penetration and strength.

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
Due to its thin, flat, knife-like shape, and lateral, bottom cutting edges, passing through tissue is very smooth thus minimizing trauma and avoiding accidental penetration of deeper tissue. Recommended for Ophthalmic procedures. (For example Cataract and Front Segment Surgery).

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
Round bodied needle with a special precision cutting point, allows easier penetration through calcified tissue. The needle is manufactured from a special stainless steel alloy that provides excellent tissue penetration and high resistance to bending. Ideal for placing heart valves in Cardiovascular Surgery.

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
Cylindrical bodied needle with a rounded point (Rome). Particularly designed for suturing soft tissue such as in liver and kidneys. This needle has a special coating that eases penetration through tissue. Due to its tip profile, it eliminates accidental glove puncture; while remaining sharp enough to penetrate fascia and muscle.

Note: This needle should not be used for skin.

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
GEOMETRY, A New Generation of Needles.

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
This new reverse cutting needle is doubly sharpened with a triangular cutting pattern. The improved geometry of the needle reduces the cutting edges angle, improving penetration and control through the tissue. Frequently used in reconstructive Plastic Surgery.

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
This new conventional needle is doubly sharpened with a triangular cutting pattern. The improved geometry of the needle reduces the angle of the cutting edges, improving the penetration and control through the tissue. Frequently used in reconstructive Plastic Surgery.

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
A taper point needle which has a cylindrical body starting at the point of assembly. The middle portion of the needle presents a special square section; designed to give superior control when using with needle holders as well as enhanced physical strength. From the proximal portion to the edge it is characterised by a conic shape (with a gradual slimming) until the needle point. The "SB" needle has a special coating that favours penetration through tissue.
Uses: General Surgery, Gastroenterology, Gynaecology and Obstetrics, Cardiovascular and Neurosurgery.

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
Cylindrical bodied needle with high precision cutting edge. The middle portion of the needle presents a special square section; it is designed to provide superior control when used in conjunction with needle holders as well as a physically superior strength. The proximal portion to the needle point has three triangular cutting edges, specifically designed for piercing through calcified tissue and atheromatous plaques with great ease. The "YB" needle has a special coating that enhances penetration through tissue.
Uses: Cardiovascular and Micro-vascular surgery.

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
A taper point needle which has a cylindrical body starting at the point of assembly. The middle portion of the needle presents a special square section; designed to give superior control when it is used with needle holders as well as enhanced physical strength. From the proximal portion to the edge, it is characterized by a conic shape (with a gradual slimming) until the needle point. The "SV" needle has a special coating that favors penetration through the tissue. The dark surface of the needle allows its easy visualization in the surgical theatre by avoiding any reflection. Very suitable in micro-vascular surgery.
Uses: Cardiovascular and Micro-vascular surgery.

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
Cylindrical bodied needle with a high precision cutting edge. The middle portion of the needle presents a square section; designed to provide superior control when used in conjunction with needle holders as well as a physically superior strength. The proximal portion to the needle point has three triangular cutting edges, specifically designed for piercing through calcified tissue and atheromatous plaques with great ease. The "YV" needle has a special coating that enhances penetration through the tissue. The dark surface of the needle allows improved visualization in the surgical theatre by avoiding any reflection. Very suitable for micro-vascular surgery.
Uses: Cardiovascular Surgery.

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
- ✦ Synthetic, absorbable, multifilament braided suture, in violet or undyed color.
- ✦ Made of a polyglycolic acid with poly (glycolide-co-lactide) (30/70) and calcium stearate coating (Violet).
- ✦ Made of a polyglycolic acid with ε-caprolactone and calcium stearate coating (Undyed).
- ✦ The tissue reactivity in microscopic form is minimal.
- ✦ Absorption occurs through progressive hydrolytic action, completed between 60 and 90 days.
- ✦ The material retains approximately 75% of its initial tensile strength by the end of two weeks, and 50% by the third week.
- ✦ **Color code:** Violet
- ✦ Used in soft tissues coaptation, ligatures, and a wide variety of surgical specialties.
See uses.

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


INITIAL STRENGTH TENSION

GENERAL CHARACTERISTICS:

The **ATRAMAT® Polyglycolic Acid** is coated with Poly(glicolide-co-lactide) (30/70) and Calcium Stearate (Violet) or with Polycaprolactone and Calcium Stearate (Undyed) in the following approximate percentages:

Polyglycolic Acid 99%
Coating 1%

CHARACTERISTICS:

- High tensile strength.
- Braided structure.
- Absorption through hydrolysis.
- Cylindrical coated multifilament thread.
- Excellent knot coupling.
- Gauge within U.S.P./E.P. guidelines.

BENEFITS:

- Excellent tissue support.
- Excellent handling.
- Predictable absorption rate and minimal tissue reaction.
- Smooth passage through the tissue.
- Providing excellent security.
- Less suture material implanted in the tissue.

When implanted in the tissue, the **ATRAMAT® Polyglycolic Acid** has minimal tissue reaction on a microscopic basis. This is comparable to other absorbable synthetic sutures of Polyglycolic Acid. However, **ATRAMAT® Polyglycolic Acid**, offers a predictable absorption rate and high tensile strength over and above the USP standard, thus providing vital artificial support to the tissue during the critical healing period. The absorption occurs through progressive hydrolytic action completed in approximately 60 to 90 days. The material retains approximately 75% of its initial tensile strength at the end of the second week, and 50% by the third week.

RECOMMENDATIONS:

As with most of the braided suture materials, it is recommended to tie a square knot with 2 or 3 additional throws. The ends of the knot should be 2 to 3 mm long depending on the surgical circumstances and the surgeon's experience.

USES:

It is indicated in all surgical procedures, where absorbable sutures are recommended. These include: General, Gastroenterology, Gynaecology, Obstetrics, Ophthalmic, Urology, Plastic surgery, and Orthopaedic surgeries.

Caution must be taken when used in elderly, malnourished or immunologically deficient patients, in which case the critical cicatrization period of the wound may be delayed.

CONTRAINDICATIONS:

Since it is an absorbable suture, it should not be used where extended tissue support is needed for an indefinite period of time.

RAW MATERIAL:

Coated Homopolymer of Glycolic Acid.

DYE:

Violet No. 2 FDC, Beige-undyed.

STERILISATION METHOD:

Ethylene Oxide.

COATING:

Poly(glycolide-co-lactide) (30/70) and Calcium Stearate (Violet).
Polycaprolactone and Calcium Stearate (Undyed).

PRESENTATION:

Box with 12, 24 or 36 sutures with one or two needles, gauges 2 to 8-0 (U.S.P.).
Box with 12, 24 or 36 pouches with pre-cut sutures, without needle, gauges 2 to 8-0 (U.S.P.).
Reels in gauges 2 to 5-0 (U.S.P.).

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- ✦ Synthetic, absorbable, multifilament braided suture, in an undyed presentation.
- ✦ Made of polyglycolic acid with polycaprolactone and calcium stearate coating.
- ✦ The tissue reactivity in microscopic form is minimal.
- ✦ Absorption occurs through progressive hydrolytic action, completed by 42 days.
- ✦ The material retains approximately 50% of its tensile strength by the end of the first week.
- ✦ **Color code:** Red.
- ✦ Frequently used in tissue that requires a short healing period.
See uses.

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Quick Absorption Polyglycolic Acid

GENERAL CHARACTERISTICS:

The **ATRAMAT® Quick Absorption Polyglycolic Acid** is coated with Polycaprolactone and Calcium Stearate in the following approximate percentages:

Polyglycolic Acid 99%
Coating 1%


CHARACTERISTICS:

- ✦ High tensile strength.
- ✦ Braided structure.
- ✦ Absorption through hydrolysis.
- ✦ Cylindrical coated multifilament thread.
- ✦ Excellent knot coupling.
- ✦ Gauge conforms to U.S.P./E.P. guidelines.

BENEFITS:

- ✦ Excellent tissue support.
- ✦ Excellent handling.
- ✦ Predictable absorption rate and minimal tissue reaction.
- ✦ Smooth passage through the tissue.
- ✦ Providing excellent security.
- ✦ Less suture material implanted in the tissue.

When implanted in the tissue, the **ATRAMAT® Quick Absorption Polyglycolic Acid** has minimal tissue reaction on a microscopic basis. This is comparable to other absorbable synthetic sutures of Polyglycolic Acid. However, **ATRAMAT® Quick Absorption Polyglycolic Acid**, offers a predictable absorption rate and high tensile strength over and above the USP standard, thus providing vital artificial support to the tissue during the critical healing period. The absorption occurs through progressive hydrolytic action completed in approximately 42 days. The material retains approximately 50% of its initial tensile strength at the end of the first week.

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
Caution must be exercised when used in elderly, malnourished or immunologically deficient patients, in which the critical cicatrization period of the wound may be delayed.

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Box with 12, 24, 36 pouches with pre-cut sutures, without needle, gauges 2 to 8-0 (U.S.P.).

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

- Synthetic, absorbable, monofilament suture, presented in a violet color or undyed.
- P-dioxanone polymer.
- Offers minimal tissue reactivity.
- Absorption occurs through slow hydrolytic action completed in approximately 180-210 days.
- The material retains approximately 70% of its tensile strength by the 4th week, and 50% by the 6th week.
- **Color code:** Silver.
- Frequently used in tissue coaptation requiring a longer healing period.
See uses.

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
SUTURE POLYDIOXANONE TENSILE STRENGTH

CHARACTERISTICS:

- High tensile strength.
- Absorption through hydrolysis.
- Monofilament with uniform diameter.
- Excellent knot coupling.
- Violet color.
- Molecular structure with soft segments.

BENEFITS:

- Excellent tissue support.
- Predictable absorption rate with minimal tissue reaction.
- Smooth passage through tissue.
- Providing excellent security.
- Excellent visibility in the operating theatre.
- Excellent handling.

INDICATIONS:

Recommended for any surgical procedures requiring a single-thread, synthetic, absorbable suture.

KNOTTING: When used with thick gauges (2, 1, 0, 2-0) it is preferable to tie both a surgeon knot and simple knot. This will allow the surgeon to adjust the tension accordingly. Once adjusted, an additional square knot is recommended according to the surgical circumstances and the surgeon experience.

When applied in thin gauges (3-0, 4-0, 5-0, 6-0, 7-0), initially two square knots are recommended with 2 or 3 additional bowknots, leaving the ends of the knots 4 to 5 mm long subject to the surgical circumstances and the surgeon experience.

USES:

General, Gastroenterology, Urology, Plastic, Gynaecology and Cardiovascular surgery.

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N/A.

PRESENTATION:

Box with 12, 24 or 36 sutures with one or two needles, gauges 2 to 7-0 (U.S.P.).

Box with 12, 24, 36 pouches with pre-cut sutures without needle, gauges 2 to 7-0 (U.S.P.).

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- ✦ Synthetic, absorbable, multifilament braided suture, in violet color or undyed.
- ✦ Made of a copolymer of glycolide acid and L-lactide [poly (glycolide-co-L-lactide)].
- ✦ The tissue reactivity in microscopic form is minimal.
- ✦ Absorption occurs through progressive hydrolytic action, completed between 56 and 70 days.
- ✦ The material retains approximately 75% of its initial tensile strength by the end of the second week, and 40% to 50% by the third week.
- ✦ **Color code:** Violet
- ✦ Frequently used in soft tissue coaptation procedures.
See uses.

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



INITIAL STRENGTH TENSION

GENERAL CHARACTERISTICS:

The **ATRAMAT® PGLA90** has the following approximate percentages:

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CHARACTERISTICS:

- High tensile strength.
- Braided structure.
- Absorption through hydrolysis.
- Cylindrical coated multifilament thread.
- Excellent knot coupling.
- Gauge within U.S.P./E.P. guidelines.
- Hermetically sealed packaging.

BENEFITS:

- Excellent tissue support.
- Excellent handling.
- Predictable absorption rate and minimal tissue reaction.
- Smooth passage through the tissue.
- Plastic coupling of the knot providing excellent security.
- Higher tensile strength with less diameter.
- Ensure hermetically sealed packaging, therefore guaranteed tightness and product sterility.

When implanted in the tissue, the **ATRAMAT® PGLA90** has minimal tissue reaction on a microscopic basis. However, **ATRAMAT® PGLA90**, offers a predictable absorption rate and high tensile strength above the USP standard, thus providing vital artificial support to the tissue during the critical healing period. The absorption occurs through progressive hydrolytic action completed in approximately 56 to 70 days. The material retains approximately 75% of its initial tensile strength at the end of the second week, and 40% to 50% by the third week.

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USES:

It is indicated in all surgical procedures, where absorbable sutures are recommended. These include: General, Gastroenterology, Gynaecology, Obstetrics, Ophthalmic, Urology, Plastic, and Orthopaedic surgery.

Note: Caution must be taken when used in elderly, malnourished or immunologically deficient patients, in which the critical cicatrization period of the wound may be delayed.

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Ethylene Oxide.

COATING:

Poly (glycolide-co-L-lactide) and Calcium Stearate.

PRESENTATION:

Box with 12, 24 to 36 sutures with one or two needles, gauges 2 to 8-0 (U.S.P.).
 Box with 12, 24 to 36 pouches with pre-cut sutures without needle, gauges 2 to 8-0 (U.S.P.).
 Reels in gauges 2 to 5-0 (U.S.P.).

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- ❖ Synthetic, absorbable, multifilament braided suture, in an undyed presentation.
 - ❖ Made of a copolymer of glycolide acid and L-lactide [poly (glycolide-co-L-lactide)].
 - ❖ The tissue reactivity in microscopic form is minimal.
 - ❖ Absorption occurs through progressive hydrolytic action, completed by 42 days.
 - ❖ The material retains approximately 50% of its tensile strength by the end of the first week.
 - ❖ **Color code:** Pink.
 - ❖ Frequently used in approaching soft tissues, skin (intradermic) and mucuose, etc., that require a short healing period.
- See uses.**

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Quick Absorption Polyglactin 910

GENERAL CHARACTERISTICS:

The **ATRAMAT PGLA90 Rapid** is composed by homopolymeric of glycolic acid with low molecular weight, braided, sterile and absorption by hydrolysis, coated with Polycaprolactone and Calcium Stearate in the following approximate percentages:

Polyglycolic Acid 90%
 L-lactide 10%
 Coating ≤1%
 [Poly (glycolide-co-L-lactide) and calcium stearate]


CHARACTERISTICS:

- ✦ High tensile strength.
- ✦ Braided structure.
- ✦ Absorption through hydrolysis.
- ✦ Cylindrical coated multifilament thread.
- ✦ Excellent knot coupling.
- ✦ Gauge conforms to U.S.P./E.P. guidelines.
- ✦ Hermetically sealed packaging.

BENEFITS:

- ✦ Excellent tissue support.
- ✦ Excellent handling.
- ✦ Predictable absorption rate and minimal tissue reaction.
- ✦ Smooth passage through the tissue.
- ✦ Providing excellent security.
- ✦ Less suture material implanted in the tissue.
- ✦ Ensure hermetically sealed packaging and therefore the stability and sterility of the product.

ABSORPTION RATE

The **PGLA90 Rapid ATRAMAT**, **Quick Absorption Polyglactin 910**, offers a predictable absorption rate and high tensile strength above the USP standard, thus providing vital artificial support to the tissue during the critical healing period. The absorption occurs through progressive hydrolytic action completed in approximately 42 days. The material retains approximately 50% of its initial tensile strength at the end of the first week.

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USES:

General, Gastrosurgery, Gynaecology and Obstetrics, Urology, Ophthalmic, Plastic, and Orthopaedic surgery.


Note: Caution must be exercised when used in elderly, malnourished or immunologically deficient patients, in which the critical healing period of the wound may be delayed.

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SUTURE PGLA90 RAPID TENSILE STRENGTH

RAW MATERIAL:

Coated Homopolymer of Glycolic Acid with low molecular weight.

DYE:

Undyed (Beige).

STERILISATION METHOD:

Gamma Radiation.

COATING:

Polycaprolactone and Calcium Stearate.

PRESENTATION:

Box with 12, 24 or 36 sutures with one or two needles, gauges 2 to 8-0 (U.S.P.).

Box with 12, 24, 36 pouches with pre-cut sutures without needle, gauges 2 to 8-0 (U.S.P.).

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- ✦ Synthetic, absorbable, monofilament suture, presented in violet color or undyed.
- ✦ Made of a copolymer of glycolide and caprolactone poly (glycolide-co-Epsilon-caprolactone).
- ✦ The tissue reactivity in microscopic form is minimal.
- ✦ Absorption occurs through progressive hydrolytic action; completed between 90 and 120 days.
- ✦ The material retains approximately 60% of its initial tensile strength by the first week.
- ✦ **Color code:** Orange soft.
- ✦ Frequently used in soft tissue coaptation.
See uses.

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

Poly (glycolide-co-Epsilon-caprolactone)

ABSORPTION RATE

SUTURE PGC25 TENSILE STRENGTH

GENERAL CHARACTERISTICS:

The **ATRAMAT® PGC25** is a copolymer of glycolide and Poly (glycolide-co-Epsilon-caprolactone) in the following approximate percentages:

Polyglycolic Acid75%
Polycaprolactone25%

CHARACTERISTICS:

- ✦ High tensile strength.
- ✦ Absorption through hydrolysis.
- ✦ Monofilament with uniform diameter.
- ✦ Excellent knot coupling.
- ✦ Violet color.
- ✦ Hermetically sealed packaging.

BENEFITS:

- ✦ Excellent tissue support.
- ✦ Predictable absorption rate with minimal tissue reaction.
- ✦ Smooth passage through tissue.
- ✦ Providing excellent security.
- ✦ Excellent visibility in the operating theatre.
- ✦ Guaranteed product sterility.

INDICATIONS:

Indicated for general approximation of soft tissue.

KNOTTING: When used with thick gauges (1, 0, 2-0) it is preferable to tie both a surgeon knot and simple knot. This will allow the surgeon to adjust the tension accordingly. An additional square knot is recommended according to the surgical circumstances and the surgeon experience.

When applied in thin gauges (3-0, 4-0, 5-0, 6-0), initially two square knots are recommended with 2 or 3 additional throws, leaving the ends of the knots 4 to 5 mm long subject to the surgical circumstances and the surgeon experience.

Damaging the suture thread with forceps, needle holders and so on, should be avoided.

USES:

General, Urology, Plastic and Gynaecological surgery.

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CONTRAINDICATIONS:

As an absorbable suture, it should not be used when the tissue support is needed for an indefinite period of time.

Note: In light of the characteristics, benefits and limitations presented above, the end user is in a position to call on his/her experience and adapt the use of the suture subject to individual surgical circumstances.

RAW MATERIAL:

Poly (glycolide-co-Epsilon-caprolactone).

DYE:

Violet No. 2 D&C, Beige-undyed.

STERILISATION METHOD:

Ethylene Oxide.

COATING:

Without coating.

PRESENTATION:

Box with 12, 24 to 36 sutures with one or two needles, gauges 2 to 6-0 (U.S.P.).
Box with 12, 24 to 36 pouches with pre-cut sutures, without needle, gauges 2 to 6-0 (U.S.P.).

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

- ✦ Animal origin, absorbable suture with twisted multifilament, in brown color.
 - ✦ Obtained from the small intestine serous layer of a healthy bovine (free of BSE and apthose fever).
 - ✦ Because it is an animal origin material, the tissue reactivity is relatively moderated.
 - ✦ Absorbed by fagositosis in approximately 60-90 days.
 - ✦ The thread keeps its tensil strength between 14 and 21 days. Specific patient conditions may cause the absorption period variation.
 - ✦ **Color code:** Brown
 - ✦ Frequently used in easy healing tissues and that do not require permanent artificial support.
- See uses.**

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CHROMIC GUT SUTURE TENSILE STRENGTH STUDY

**PRESERVATION SOLUTION**

Isopropyl alcohol	89%
Distilled water	9%
Triethanolamine	2%

Total 100%

CHARACTERISTICS:

- High purity collagen between 97% and 98%.
- Chromicizing process before twisting.
- Uniform calibration and polishing, monofilament appearance.
- Sterilized by gamma rays of Cobalt 60

BENEFITS:

- Minimum tissue reactivity.
- Uniform chroming along the thread length, providing a reliable tissue support during the absorption.
- Provides smooth passage through the tissue, minimizing tissue drag.
- Complete and reliable sterilization, preserving the collagen natural characteristics.

USES:

General, Gynaecology, Obstetrics, Ophthalmic, Urology, and other surgery.

Note: Caution must be taken when used in elderly, malnourished or immunologically deficient patients, in which the critical healing period of the wound may be delayed. In light of this, the surgeon may exercise with discretion.

CONTRAINDICATIONS:

As a natural origin, absorbable suture, it should not be used where tissue support for an extended period of time is required.

For further technical information please refer to the instructive attached in each product box.

CONSIDERATIONS:**WHAT TO DO.**

HANDLE the surgical gut only when wet.

STRETCH the surgical gut immediately after taking it out of the envelope pulling should be uniform, firm and only one soft pull.

STRAIGHT the suture and the needle, holding the thread approximately 2 inches (5 cm) from the joint and the distal end of the needle.

HUMIDIFY the surgical gut when the material has dried. Humidifying it only in a saline solution or in sterile water for not more than 4 to 6 seconds, just before giving it to the surgeon.

HANDLE carefully the surgical gut only for the necessary time.

WHAT NOT TO DO.

DO NOT handle the surgical gut when dry.

DO NOT try to stretch the surgical gut when it is already dry.

DO NOT expose the surgical gut to unnecessary stress, twisting or hard pulling of the thread.

DO NOT pass the gloved fingers on the surgical gut thread while stretching it to remove the memory.

DO NOT pull the suture-needle joint.

DO NOT dip the surgical gut in solution for more than 6 seconds.

DO NOT press the surgical gut with instruments such as forceps or needle supports, etc.

RAW MATERIAL:

Serous layer of the small intestine from healthy bovine cattle (free of BSE and aphtose fever).

CHROMATIZED:

Chemical bath with chromic salts solutions.

STERILISATION METHOD:

Gamma Radiation Cobalt 60.

COATING:

None.

PRESENTATION:

Box with 12, 24 or 36 sutures with one or two needles, gauges 2 to 6-0 (U.S.P.).

Box with 12, 24 or 36 pouches with pre-cut sutures without needle, gauges 2 to 6-0 (U.S.P.).

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- ✦ Animal origin, absorbable suture with twisted multifilament, in yellow color.
 - ✦ Obtained from the small intestine serous layer of a healthy bovine (free of BSE and apthose fever).
 - ✦ Absorbed by fagositosis in approximately 50-70 days.
 - ✦ The thread keeps its tensil strength between 7 and 14 days. Specific patient conditions may cause the absorption period variation.
 - ✦ **Color code:** Yellow
 - ✦ Frequently used in easy healing tissues and that do not require permanent artificial support.
- See uses.**

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**CHARACTERISTICS:**

- High purity collagen between 97% and 98%.
- Uniform calibration and polishing, monofilament appearance.
- Sterilized by gamma rays of Cobalt 60.

BENEFITS:

- Minimum tissue reactivity.
- Provides smooth passage through the tissue, minimizing tissue drag.
- Complete and reliable sterilization, preserving the natural characteristics of collagen.

USES:

General Surgery, Gynaecology, Obstetrics, Ophthalmic, etc.

Note: Caution must be taken when used in the elderly, malnourished or immunologically deficient patients, in which the critical healing period of the wound may be delayed. In light of this, the surgeon may exercise with discretion.

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For more technical information please refer to the instructive attached in each product box.

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
DO NOT expose the surgical gut to an unnecessary stress, twisting or hard pulling of the thread.

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ABSORPTION RATE**PLAIN GUT SUTURE TENSILE STRENGTH STUDY (IN VIVO)****PRESERVATION SOLUTION**

Isopropyl alcohol 89%
Distilled water 9%
Triethanolamine 2%

Total 100%

RAW MATERIAL:

Serous layer of the small intestine from healthy bovine cattle (free of BSE and apthose fever).

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- Monofilament, synthetic, non-absorbable suture, presented in black, blue or undyed.
- Obtained from extrusion of polyamide (6.0 and 6.6) with a uniform cylindrical diameter.
- Tissue reaction is minimal.
- Nylon is a non-absorbable material that with time becomes encapsulated by connective tissue.
- It is known that the thread mass diminishes, approximately by 10% a year due to the rupture of chemical links in the thread as a result of hydrolytic action.
- **Color code:** Green
- Normally used when confronting tissue in Neurological, Ophthalmic and Plastic Surgery. **See uses.**

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GENERAL CHARACTERISTICS:

Its tensile strength and controlled elongation ensures the approximation of tissue borders. In addition, the sutures smooth passage through the tissue guarantees perfect healing. These physical attributes coupled with minimal tissue reaction, constitute an ideal suture for confronting a variety of tissue. Our **ATRAMAT® Nylon** materials are accompanied by wide range of non-traumatic needles, precisely sharpened, allowing for favorable cosmetic results. (See the needle chart)

CHARACTERISTICS:

- ✦ Synthetic origin (prime quality raw material). Inert.
- ✦ Excellent knot coupling.
- ✦ Monofilament with uniform gauge.
- ✦ Hermetically sealed packaging.
- ✦ Black and blue color.

BENEFITS:

- ✦ Low tissue reaction (upon implantation).
- ✦ Providing excellent knot security.
- ✦ Smooth passage through tissue therefore, avoiding tissue drag.
- ✦ Guaranteed product sterility.
- ✦ Excellent visibility in surgical field.

USES:

NYLON (Blue, Black)

Plastic, Gynaecology, Obstetrics, Orthopaedic and General surgery.

NYLON (Black)

Ophthalmology, Neurology, Micro-vascular surgery.

NYLON (Undyed)

Plastic and Reconstructive surgery.

RAW MATERIAL:

Polyamide 6.0 or 6.6.

DYE:

Black Hematine, Blue No. 2 FDC.

STERILISATION METHOD:

Gamma Radiation or Ethylene Oxide.

COATING:

None.

PRESENTATION:

Box with 12, 24 to 36 sutures with one or two needles, gauges 2 to 11-0 (U.S.P.).

Box with 12, 24 to 36 pouches with precut sutures without needle, gauges 2 to 11-0 (U.S.P.).

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- ❖ Synthetic, non-absorbable, monofilament suture.
- ❖ Presented in a blue color.
- ❖ Extruded as a filament with a computer controlled diameter.
- ❖ The tissue reaction is minimal.
- ❖ The Polypropylene is extraordinarily stable; ideal for fulfilling its purpose for permanent tissue support without compromising its tensile strength.
- ❖ **Color code:** Intense blue.
- ❖ Frequently used to confront tissue in specialised areas. Plastic and Cardiovascular procedures being among the most important.
See uses.

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**GENERAL CHARACTERISTICS:**

Its high tensile strength and controlled elongation ensures the approximation of tissue borders.

The smooth passage through the tissue ensures perfect cicatrisation. In addition, its low tissue reaction makes it ideal for cuticular use, and since it is able to withstand repeated bending for long periods of time, it is also ideal in Cardiovascular surgery. Its physical and biological attributes offer the following:

CHARACTERISTICS:

- ❖ Synthetic origin Inert. Prime quality raw material.
- ❖ Excellent plastodeformation.
- ❖ Monofilament (uniform gauge).
- ❖ Hermetically sealed packaging.
- ❖ Non-absorbable.
- ❖ Resistance to repeated bending (flexion).

BENEFITS:

- ❖ Low tissue reaction, ideal biocompatibility.
- ❖ Adequate coupling of the knot providing excellent knot security.
- ❖ Smooth passage through tissue, thus eliminating tissue drag.
- ❖ Guaranteed product sterility.
- ❖ Comprehensive and extended tissue support.
- ❖ Ideal for Cardiovascular procedures requiring support for the heart's propulsion, either in heart prosthesis or vascular anastomosis.

EASY PENETRATION WITHOUT FILTRATION (Cardiovascular Surgery)

Our suture material has non-traumatic needles, precisely sharpened and specially designed to provide favourable cosmetic results. Furthermore, our needles have special bending (flexion) properties, resistance and hardness.

The Rockwell hardness of the needle assembled to Polypropylene is higher than the conventional needles. This allows excellent penetration and passage through calcified tissue or atheromatous plaques.

ATRAMAT® have special needles for cardiovascular use, which provide easy penetration of the tissue, keeping a uniform diameter in the thread-to-needle ratio. This helps to avoid filtration, and is ideal for vascular use.

USES:

Cardiovascular, Plastic and Reconstructive, Gynaecology, Obstetrics, Orthopaedic and General surgery.

RAW MATERIAL:

It is an isotactic estereoisomer (fiber forming agent) of Polypropylene Homopolymer.

DYE:

Blue No. 2 FDC.

STERILISATION METHOD:

Ethylene Oxide.

COATING:

Without coating

PRESENTATION:

Box with 12, 24 to 36 sutures, individual presentation pledgets, in gauges 3-0 and 2-0 (U.S.P.).

Box with 12, 24 to 36 sutures with one or two needles, in gauges from 2 to 10-0 (U.S.P.).

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

POLYPROPYLENE

Ideal for Cardiovascular procedures

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- ❖ Synthetic, non-absorbable, multifilament braided suture.
- ❖ Presented in a green or white color.
- ❖ Polyester made of Terephthalate with or without coating.
- ❖ Its non-absorbable, synthetic origin means it has minimal tissue reaction.
- ❖ Used in tissue coaptation due to its characteristically high tensile strength.
- ❖ **Color code:** Orange
- ❖ Frequently used in speciality surgery, including Cardiovascular and Ophthalmic procedures.
See uses.

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**CHARACTERISTICS:**

- ✦ Synthetic origin. Inert.
- ✦ Resistant to repeated flexion.
- ✦ Braided multifilament.
- ✦ Non-absorbable.
- ✦ Green or white color.
- ✦ Excellent knot coupling.
- ✦ Hermetically sealed packaging.

BENEFITS:

- ✦ Minimal tissue reaction.
- ✦ Ideal for use in Cardiovascular Surgery. (Heart Prosthesis)
- ✦ Superior tensile strength.
- ✦ Extended tissue support.
- ✦ Good visibility in the surgical theatre.
- ✦ Providing excellent knot security.
- ✦ Guaranteed product sterility.

The Rockwell hardness of the needle assembled to the Polyester suture has a higher quality than the one used in conventional needles. This provides excellent penetration; easing the way through calcified tissue or atheromatous plaques.

USES:

It is indicated in procedures where the use of a synthetic, braided, non-absorbable suture is recommended. These areas include: Cardiovascular, Orthopaedic, Ophthalmic and General surgery.

RAW MATERIAL:

Polyester fibres prepared from Terephthalic acid or its esters and ethylene glycol.

DYE:

Green No. 6 FDC, White.

STERILISATION METHOD:

Gamma Radiation or Ethylene Oxide.

COATING:

Silicone or teflon.

PRESENTATION:

Box with 12, 24 to 36 sutures with one or two needles, with or without coating, gauges 5 to 7-0 (U.S.P.).

Pledgets, presentation in 10 suture kit, thread colors: green or white, gauges 2-0 to 3-0 (U.S.P.).

Pledgets, individual presentation, gauges 0 to 5-0 (U.S.P.), thread colours: green or white.

Box with 12, 24 to 36 pouches with pre-cut sutures without needle, with more than two threads, gauges 2-0 to 3-0 (U.S.P.).

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- ❖ Synthetic, non-absorbable, monofilament suture, in blue or yellow color.
 - ❖ The PVDF is a Polyvinylidene Fluoride polymer.
 - ❖ Extruded as a filament with a computer controlled diameter.
 - ❖ The PVDF in vivo is extraordinarily stable, ideal for fulfilling its purpose for permanent tissue support; without compromising its tensile strength.
 - ❖ The chemical structure of PVDF, provides greater resistance to repeated bending.
 - ❖ The tissue reaction is minimal.
 - ❖ **Color code:** Blue
 - ❖ Frequently used to confront tissue in specialized areas: Cuticular and Cardiovascular procedures are among the most important.
- See uses.**

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GENERAL CHARACTERISTICS:

The surgical suture monofilament **PVDF ATRAMAT®** has a chemical structure that provides greater resistance to repeated bending and stability compared with polypropylene suture and also allow a knot coupling without breakdown, giving a good knot security, if the applicable knotted technique for monofilament suture is used.

Its high tensile strength and controlled elongation ensures the approximation of tissue borders until healing.

Due to its smooth passage through the tissue ensures perfect healing. In addition, its low tissue reaction makes it ideal for cuticular use and since it resists repeated bending for long periods of time, the **PVDF ATRAMAT®** is a biocompatible and anti-thrombogenic material ideal for Cardiovascular Surgery. Due to all the physical and biological characteristics it offers:

CHARACTERISTICS:

- ❖ Synthetic origin inert. Prime quality raw material.
- ❖ Excellent knot coupling.
- ❖ Monofilament with millimetrical uniform gauge.
- ❖ Blue and yellow color.
- ❖ Adequate ratio between needle and thread diameter.
- ❖ Non-absorbable.
- ❖ Resistance to repeated bending.
- ❖ Needle protection support.
- ❖ Precise sharpness premium needles.
- ❖ Hermetically sealed packaging.

BENEFITS:

- ❖ Low tissue reaction, ideal biocompatibility and anti-thrombosis.
- ❖ Adequate coupling of the knot providing excellent knot security.
- ❖ Smooth passage through the tissue, avoiding tissue drag.
- ❖ Excellent visibility in surgical field as well as versatile applications.
- ❖ Ideal for cardiovascular procedures, avoid filtration.
- ❖ Adequate and extended tissue support.
- ❖ Comprehensive and extended tissue support. Ideal for cardiovascular procedures, supporting the heart propulsion, either in heart prosthesis or vascular anastomosis avoid filtration.
- ❖ Precise sharp needles protection.
- ❖ Excellent penetration with minimal trauma and cosmetic results.
- ❖ Ensures hermetically sealed packaging, therefore guaranteed product sterility.

EASY PENETRATION WITHOUT FILTRATION

Our suture material has non-traumatic needles, precisely sharpened and specially designed to provide favorable cosmetic results. Furthermore, our needles have special bending properties, resistance and hardness.

The Rockwell hardness of the needle assembled to PVDF is higher than the conventional needles. This allows excellent penetration and passage through calcified tissue or atheromatous plaques.

We have special needles for cardiovascular use, which provide easy penetration of the tissue, keeping a uniform diameter in the thread-to-needle ratio. This helps to avoid filtration, and is ideal for vascular use.

USES:

Cardiovascular, Plastic Reconstructive, Cuticular, Gynaecology, Obstetrics, Orthopaedics and General surgery.

RAW MATERIAL:

Polyvinylidene Fluoride Polymer.

DYE:

Blue and Yellow.

STERILISATION METHOD:

Ethylene Oxide.

COATING:

None.

PRESENTATION:

Box with 12, 24 and 36 sutures with Pledgets, individual presentation, gauges 4-0 to 2-0 (U.S.P.).
Box with 12, 24 and 36 pouches with pre-cut sutures with one or two needles, gauges from 2 to 8-0 (U.S.P.).

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CHARACTERISTICS:

- ✦ Prime quality raw material.
- ✦ Threaded multifilament.
- ✦ Non-absorbable, natural origin.
- ✦ Silicone Coating.
- ✦ Blue, white color (virgin silk), with / without silicone coating.
- ✦ Hermetically sealed packaging.

BENEFITS:

- ✦ Prime quality guaranteed silk.
- ✦ The high technology used in the threading guarantees a compact thread with a higher tensile strength.
- ✦ Longer tissue support.
- ✦ Smooth passage through tissue.
- ✦ Versatile applications.
- ✦ Guaranteed product sterility.

INDICATIONS:

Indicated in all surgical procedures provided the limitations of this material are known.

USES:

General, Gastroenterology, Ophthalmology, Gynaecology, and Obstetric surgery.

CONTRAINDICATIONS:

Since it is a non-absorbable suture, it is NOT recommended in Urological procedures as nor biliary tracts as this may lead to the formation of stones around the suture.

RAW MATERIAL:

It is made with natural silk threads produced by the Bombyx Mori I Silk Worm.

DYE:

Black Hematine Silk, Blue-blue Silk No. 2 FDC, Virgin Silk-without color.

STERILISATION METHOD:

Ethylene Oxide or Gamma Radiation.

COATING:

Silicone.

PRESENTATION:

Box with 12, 24 to 36 sutures with one or two needles, gauges 4 to 8-0 (U.S.P.).
Box with 12, 24 to 36 pouches with pre-cut sutures without needle, gauges 4 to 6-0 (U.S.P.).
Reels in gauges 4 to 5-0 (U.S.P.).

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


Atramat

Silk

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Excellent Visibility in Surgical Field

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
- ✦ Monofilament, stainless steel, non-absorbable suture.
- ✦ Presented in a natural color.
- ✦ The tissue reaction is minimal.
- ✦ Permanent tensile strength.
- ✦ **Color code:** Yellow.
- ✦ Frequently used in Orthopaedic Surgery including Sternum closure.
See uses.

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Stainless Steel

CHARACTERISTICS:

- ✦ Stainless steel.
- ✦ Non-absorbable.
- ✦ Monofilament.
- ✦ Permanent high tensile strength.
- ✦ Ductile and malleable.
- ✦ Special non-traumatic needles design.
- ✦ Hermetically sealed packaging.

BENEFITS:

- ✦ Inert (compatible with tissue).
- ✦ Extensive tissue support.
- ✦ Smooth passage through tissue.
- ✦ Excellent support to tissue (bone structures).
- ✦ Good handling.
- ✦ Excellent penetration in hard tissues.
- ✦ Guaranteed product sterility.

USES:

Once the user is aware of the characteristics, benefits and limitations of the stainless steel suture, and if he/she is familiar with good surgical practices, this product is adequate for procedures that require a highly resistant material for example, in Orthopaedics and closure of the sternum.

RAW MATERIAL:
Stainless Steel.

DYE:
N / A.

STERILIZATION METHOD:
Gamma Radiation or Ethylene Oxide.

COATING:
N / A.

PRESENTATION:
Box with 12 or 48 sutures with needles in gauges 2, 5 and 6 (U.S.P.).

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